

AAC app cheat sheet

Speak for Yourself

Hold That Thought: Use to say a message that you saved

Edit: Use this feature to change/add a cell

Search Icon: Search for available words and follow path to word

Action words are color coded green on the second level

Thought Bubble Icon: Use to save a message to say later

Whisper Mode: Turns off word by word speaking to allow private communication

Open/Close: Use to easily select which words to show or hide

Babble: Use to open everything to give user opportunity to explore and babble

Settings: Use this menu to change voices, select vocabulary templates, back up and more.

Lock: Use to lock editing feature. Go to iOS Settings -> SFY to unlock

Message Window: Press to say message

Arrow X: Delete one word; Circle X: Delete whole message

Sharing: Share contents of message window with connected accounts (text, email, and social media)

Sharing icon (box with up arrow) will appear when there is a message in the window

QWERTY: Access to a keyboard.

Speak for Yourself application from Speak for Yourself, LLC

Forbes AAC
 181 Illinois Ave. South
 Mansfield, OH 44905
 419.589.7688
 fax 419.589.5146



forbesaac.com