

Use Guided Access on your ProSlate

Guided Access helps you stay focused on a task by temporarily restricting the ProSlate to a single app, and allowing you to control which app features are available. You can do any of the following:

- Disable areas of the screen that aren't relevant to a task, or areas where an accidental gesture might cause a distraction
- Disable the ProSlate's hardware buttons
- Limit how long someone can use the app

Set up Guided Access

1. Go to Settings > Accessibility > Guided Access, then turn on Guided Access.
2. Adjust any of the following:
 - Passcode Settings: Tap Set Guided Access Passcode, then enter a passcode.
 - You can also turn on Face ID or Touch ID as a way to end a Guided Access session.
 - Time Limits: Play a sound or speak the time remaining before a Guided Access session ends.
 - Accessibility Shortcut: Turn the shortcut on or off during Guided Access sessions.
 - Display Auto-Lock: Set how long it takes the ProSlate to automatically lock during a Guided Access session.

Start a Guided Access session

1. Open the app you want to use.
2. Turn on Guided Access using accessibility shortcuts.
3. Circle any areas of the screen you want to disable. Drag the mask into position or use the handles to adjust its size.
4. Tap Options, then turn on or off any of the following:
 - Top Button (or Sleep/Wake Button)
 - Volume Buttons
 - Motion (to prevent ProSlate from switching from portrait to landscape or from responding to other motions)
 - Keyboards
 - Touch
 - Time Limit
5. Tap Start.

End a Guided Access session

Do any of the following:

- On the ProSlate with a Home button: Double-click the Home button, then unlock with Touch ID (if enabled). Or triple-click the Home button, then enter the Guided Access passcode.
- On an ProSlate with Face ID: Double-click the top button, then unlock with Face ID (if enabled). Or triple-click the top button, then enter the Guided Access passcode.

Forbes AAC

181 Illinois Ave. South
Mansfield, OH 44905

phone 419.589.7688

fax 419.589.5146



forbesaac.com