

AAC app cheat sheet

Grid for iPad: Super Core 50

Magic Wand: Conjugates chosen word

Jump Back: Jumps back one page

Speak: Use to speak message

... (Three Dots): Opens Edit grids, Find word, Grid set manual, and Settings

jump home

Jump Home: On other pages this will take

Purple Row: Question words

Gray Columns: Common little words

Blue Columns: Fringe words that when selected occupy the blue columns leaving core words visible

Green Columns: Action words

Blue Cells: Descriptive words

Yellow Column: Pronouns

NOTE: Cells with a + sign link to more words in that category

Spelling/ABC: Open QWERTY keyboard; Contains word search feature

Options

GRIDS

- Edit grids
- Find word
- Grid set manual

SETTINGS

- Settings

forbesaac.com/teletherapy-assessment-program

Super Core 50[®] is a trademark of Tobii Dynovox

Forbes AAC
 181 Illinois Ave. South
 Mansfield, OH 44905
 419.589.7688
 fax 419.589.5146

