

# AAC app cheat sheet

## TouchChat WordPower 42

**Vocab:** Copy and select other TouchChat page sets here

**Pronouns:** Common pronouns are coded yellow. More people words can be found in the people category

**Message Window:** Press to say message

**Menu:** Go here to change look of the board, edit cells and pages, and more.

**Arrow X:** Delete whole message

**ABC/123:** Access to a keyboard and the FIND WORD feature

**Home:** The yellow home button on other pages will always bring you back to your home page

**Fringe Words:** Words organized by category. The single and double arrows mean the cell links to another page

**Extra Words:** Connecting words, little words, and transition words

**Action Words:** Common action words coded green

**Clear:** Erases what is in the message window

Vocab						Menu
<b>PEOPLE</b> 	<b>QUESTIONS</b> ?	<b>PLACES</b> 	<b>SOCIAL</b> 	<b>TIME</b> 	<b>GROUPS</b> 	<b>ABC 123</b> 
I	can	drink	eat	feel	good	DESCRIBE 
my	do	finish	help	go	a	more 
it	have	like	listen	play	the	that 
you	is	watch	work	clear	EXTRA And It The WORDS Up	
your	don't	.	want			

[forbesaac.com/teletherapy-assessment-program](http://forbesaac.com/teletherapy-assessment-program)

TouchChat® is a trademark of PRC-Salttillo

**Forbes AAC**  
181 Illinois Ave. South  
Mansfield, OH 44905  
419.589.7688  
fax 419.589.5146



[forbesaac.com](http://forbesaac.com)