

# Sixteen plus hours of run time for an AAC device!

Listening to feedback from the field, we were tasked with solving a problem: **more run time for AAC devices**. The design engineering team at Forbes AAC has created the XTNDR™ Auxiliary Battery Pack that gives users an **amazing sixteen hours** of normal run time and more!

The XTNDR wheelchair mount battery pack provides WinSlate or ProSlate users with power and freedom without adding any bulk or weight to their AAC devices.

XTNDR is an **industry first** design concept that allows an external battery pack to be mounted directly to a wheelchair mount tube. Mount it to the horizontal crossbar behind the device or mount it to the vertical tube affixed to the chair. Or, it can be popped into a back pack for easy storage.

The battery pack features cord management retainers to take up any slack and prevent annoying snags or hang ups. XTNDR charges a WinSlate using a magnetic break-away charge cable or a ProSlate using a standard charging cable.



Use XTNDR in a horizontal or vertical configuration

Integral clamp system for 7/8" mount tube



[forbesaac.com/xtndrbatterypack](http://forbesaac.com/xtndrbatterypack)

**Forbes AAC**  
181 Illinois Ave. South  
Mansfield, OH 44905

phone 419.589.7688  
fax 419.589.5146



[forbesaac.com](http://forbesaac.com)